



Juzo Night

Nighttime maintenance for effective treatment of lymphedema



Juzo Think Green. Live Green.

We're committed to reducing our carbon footprint for a greener planet.

Wind Turbine



Solar Array



Prairie Land Grass Landscaping



Motion Sensor Lights



Hands-Free Faucets & Dryers



Recycling Centers



Juzo USA
3690 Zorn Drive
Cuyahoga Falls, OH 44223
www.juzousa.com
support@juzousa.com




Juzo Canada Ltd.
1100 Burloak Dr. Ste. 300
Burlington, ON L7L 6B2
www.juzo.ca
support@juzo.ca


©Juzo, Slippies, Slippie Gators, Keep Moving & Signature Collection are registered trademarks of Julius Zorn, Inc. • Lycra Premium Stretch Fibers is a registered trademark of the Lycra Company. Tyvek is a registered trademark of Du Pont de Nemours and Company. • X-STATIC is a registered trademark of Noble Fiber Technologies, Inc. © 2022 Julius Zorn Inc.

Juzo Night

Uniform, comfortable compression for lymphedema maintenance. Juzo Night is the perfect complement to circular and flat knit compression garments worn during the day, preventing edema from rebounding during the night.



24 Hour therapy



For rest or sleep




Comfortable, convenient




Removable fasteners

Easy to put on and take off



Thumb guard

Keeps edema from shifting




Foam channels

Follow lymphatic pathways, promoting lymphatic drainage



Anatomical angle

Prevents bunching and skin irritation, protects pressure points



Cooling material

Wicks moisture, reduces sweating





Easy to wear – anywhere

Juzo Night is patient-friendly and effective when day wear compression garments are removed for resting or sleep. Juzo Night is lightweight, comfortable and easy to put on and take off, encouraging patient self-care and compliance anywhere.

Successful self-management

- + Lining and touch fasteners make it easier to put on
- + Individually adjustable; can be adjusted at any time
- + Adjusts to fluctuations in circumference that occur over the course of treatment



Lining


Reduces slipping, keeps the garment in place



Perfect Night Garment



A good night's sleep or afternoon nap positively impacts patient well-being. Juzo Night provides ease-of-use and exceptional comfort, and reduces sensations of tension and pressure caused by edema.

Use of Juzo lymphology products during complex manual decongestive therapy



Everything going the right way

Designed to support lymphatic pathways, Juzo Night's stitched foam channels promote lymphatic drainage, ensure decongestion, and reduce swelling.



Juzo Night, A Nighttime Garment



Indications, contraindications and material composition

Scan the QR code to see usage and care instructions or go to: juzousa.com/juzonight or juzo.ca/juzonight



Color



Washable in warm water



Dryer-safe

Size chart

Juzo Night · Nighttime Garment - 6872

	Circumference c	Standard sizes · Measurements in cm		
		2 / S	3 / M	4 / L
Length L				
LCG	cG	31.0 – 35.0	34.0 – 38.0	37.0 – 41.0
	cF	29.0 – 32.0	32.0 – 35.0	35.0 – 38.0
	cE	26.0 – 29.0	29.0 – 32.0	32.0 – 35.0
	cD	23.0 – 26.0	26.0 – 29.0	29.0 – 32.0
	cC	15.0 – 17.5	17.5 – 20.0	20.0 – 22.5
Length	regular	42.0 – 47.0	42.0 – 47.0	42.0 – 47.0
	LCG long	47.0 – 51.0	47.0 – 51.0	47.0 – 51.0

• Specify right or left

PCSZ-01-72

