

# Custom Measurement Form for Stump Shrinkers

Phone: 1 800 222-4999

Fax: 1 800 645-2519



## Account Information

Account Number \_\_\_\_\_ Date \_\_\_\_\_

Account Name \_\_\_\_\_ Contact \_\_\_\_\_

Phone \_\_\_\_\_ Fax \_\_\_\_\_

Patient ID \_\_\_\_\_

Prescribing Physician \_\_\_\_\_

## Circular Knit Compression Garment Information (without seam)

Compression:  20-30 mmHg    Model:  Dynamic (Varin Soft-in)  
 30-40 mmHg     Dynamic Silver (Varin Soft-in)

Style:  C-F below knee     F-G above knee     H-D Pirogoff  
 Silicone Border     Hip Attachment

K-T Body Part (worn with F-G)     Slip-on     Hook and Loop Closure  
 20-30 mmHg     30-40 mmHg

Extremity:  Right     Left    Quantity: \_\_\_\_\_  Pieces


Comments: \_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

## Circumferences



## Lengths

- S-T .....
- S-H .....
- S-K } .....
- S-G } .....
- S-F .....
- S-E .....
- S-D .....
- S-C .....
- S-B1 .....
- S-B .....

Measure circumferences beginning 5 cm from the distal end of the stump.

**S** = Distal end of the stump. Measure all lengths from 'S'.