

Custom Measurement Form for Stump Shrinkers

Phone: 1 800 222-4999

Email: customercare@juzousa.com



Account Information

Account Number _____ Date _____

Account Name _____ Contact _____

Phone _____ Fax _____

Patient ID _____

Prescribing Physician _____

Circular Knit Compression Garment Information (without seam)

Compression: 18-21 mmHg Model: Dynamic
 23-32 mmHg Dynamic Silver

Flat Knit Compression Garment Information (with seam)

Compression: 18-21 mmHg Model: Expert
 23-32 mmHg Expert Cotton


Style: C-F below knee F-G above knee H-D Pirogoff
 Silicone Border Hip Attachment

K-T Body Part (worn with F-G) Slip-on Hook and Loop Closure
 18-21 mmHg 23-32 mmHg

Extremity: Right Left Quantity: _____ Pieces

Comments: _____

PCSZ-01-18b



Measure circumferences beginning 5 cm from the distal end of the stump.

S = Distal end of the stump. Measure all lengths from 'S'.